

GLOVE DISPENSER - plastic



Product Code	25430
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279254303
PARAF	987476619

Description

Glove dispenser - plastic

Holds one box of 100 gloves.
Made of unbreakable washable thermoplastic material.
Fixed base and removable cover once the box of gloves needs to be replaced.
Can be easily hung to the wall and into ambulance in vertical or horizontal position.
Compatible with box of gloves up to 25.5x7xh 12 cm.

Weight: 322 g
Size: 25.7x10xh 13.4 cm