

GLOVE DISPENSER - single - stainless steel



Product Code	25431
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279254310

Description

Glove dispenser - single - stainless steel

Hold 1 standard box of 100 gloves.
Made of AISI 304 stainless steel.
Autoclavable at 121 to 134 °C.
Can be easily hung to the wall and into ambulances.
Screw included.

Size: 25x9.5xh 13.5 cm