



Product Code	54315
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

CH-REP	Yes
--------	-----

RDM (NSIS)	2695840
CND	Z12050482
EAN/UPC	6971049341987
GMDN	41651

## Description

AI ECG TRACKER SOFTWARE for professional use - optional

AI-ECG tracker is an advanced professional dynamic ECG diagnostic analysis system. It is certified as a software for medical devices and has complete dynamic ECG analysis function, suitable for the evaluation of adult ECG data and detection of arrhythmia-like heart abnormalities.

A cardiologist can edit, modify and confirm the diagnosis.

Diagnosis accuracy rate of 95.2%. It supports template editing, atrial flutter and fibrillation analysis, HRV analysis, heart rate oscillation analysis, P-wave superposition analysis, ST-segment analysis and other advanced analysis functions.