

# GLOVE DISPENSER - triple - stainless steel



Product Code	25432
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279254327
---------	---------------

PARAF	987469525
-------	-----------

## Description

Glove dispenser - triple - stainless steel

Hold 3 standard boxes of 100 gloves.

Made of AISI 304 stainless steel.

Autoclavable at 121 to 134 °C.

Can be easily hung to the wall and into ambulances.

Screw included.

Suitable for S - M - L box of gloves together

Size: 25.2x9.5xh 37.7 cm