

"NOVA" LED TORCH - white



Product Code	25440
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
EAN/UPC	8023279254402
GMDN	12347

Description

Professional penlights, also suitable for diagnostic, with bulb or white LED with quick clip shifter mechanism or push button (25634). Argenta and Elegance made of anodized aluminium, Nova in topquality ABS plastic housing. Supplied with 2x1.5V AAA batteries.

Nova torch
Light: LED
Light intensity: 25,000 m
Material: plastic
Supplied with batteries: 2 x AAA
Colour: white
ON/OFF by: Pocket clip

Made in Germany.