

"ARGENTA" LED TORCH - metal - blue



Product Code	25442
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
EAN/UPC	8023279254426
GMDN	12347

Description

Professional penlights, also suitable for diagnostic, with bulb or white LED with quick clip shifter mechanism or push button (25634). Argenta and Elegance made of anodized aluminium, Nova in topquality ABS plastic housing. Supplied with 2x1.5V AAA batteries.

Argenta torch
Light: LED
Light intensity: 25,000 m
Material: metal
Supplied with batteries: 2 x AAA
Colour: blue
ON/OFF by: Pocket clip

Made in Germany.