

SILICONE TUBE 5x8 mm - 1.5 mm thick



| | |
|---------------|----------------|
| Product Code | 25460 |
| Unit of sale | roll of 30 m |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2207976 |
| CND | V9099 |
| EAN/UPC | 8023279254600 |
| PARAF | 986719021 |
| GTIN | 08023279254600 |
| GMDN | 35917 |

Description

SILICONE TUBES for external medical use - silicone 100%

Resistant to hot air up to 200°C.

Infra-red vulcanization. Hours treatment at 200°C temperature in special oven. Anti-adherent, avoid incrustation and coagulation.

Medical Grade - Transparent - Odourless

Made in European Union according to: US PHARMACOPOEIA XXI EUROPEAN PHARMACOPEA DIN 58367 - DIN 58362 - DIN 13098- FDA 177.2600 - BGA XV part A and B - BGA IX part B

Made in Spain.

Technical Specifications

Raw material: HCR

Colour: translucent

Hardness: 60° ±5% Shore - DIN 53505

Density: 1.17 g/ccm - DIN 53479A

Stretch: 490/530% - DIN 53504S1

Tear resistance: 34-35 N/mm - ASTM D 624 B