

# SILICONE TUBE 7x10 mm - 1.5 mm thick



Product Code	25462
Unit of sale	roll of 30 m
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2207976
CND	V9099
EAN/UPC	8023279254624
PARAF	987473218
GTIN	08023279254624
GMDN	35917

## Description

SILICONE TUBES for external medical use - silicone 100%

Resistant to hot air up to 200°C.

Infra-red vulcanization. Hours treatment at 200°C temperature in special oven. Anti-adherent, avoid incrustation and coagulation.

Medical Grade - Transparent - Odourless

Made in European Union according to: US PHARMACOPOEIA XXI EUROPEAN PHARMACOPEA DIN 58367 - DIN 58362 - DIN 13098- FDA 177.2600 - BGA XV part A and B - BGA IX part B

Made in Spain.

## Technical Specifications

Raw material: HCR

Colour: translucent

Hardness: 60° ±5% Shore - DIN 53505

Density: 1.17 g/ccm - DIN 53479A

Stretch: 490/530% - DIN 53504S1

Tear resistance: 34-35 N/mm - ASTM D 624 B