

S/S STRAIGHT ARTERY FORCEPS - foot print fantasy - 16 cm



Product Code	20549
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1630372
CND	L031306
EAN/UPC	8023279205497
GTIN	08023279205497
GMDN	62468

Description

Stainless steel nurse straight artery forceps - 16 cm

- Model: foot print fantasy

Powder coated instruments are not autoclavable, can be sterilized in water at maximum 110°C