

COMPATIBLE PAEDIATRIC PADS for defibrillator Medtronic Physio Control



Product Code	55004
Unit of sale	5 pairs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
RDM (NSIS)	2469016
CND	C020401
GTIN	18033695728262
GMDN	47055

Description

Pediatric Compatible pads - disposable
(CPR guidelines state 8 years as limit age for pediatric treatment, corresponding to about 25 kg according to standard growth curves)

These pads have an energy reduction system: defibrillator output of 150J (33606) or 200J (33603) is reduced to 50J.
International compliance to regulation as above adult compatible pads.
Latex free.
Made in Italy.

Multilingual envelope: GB, FR, IT, ES, PT, DE, NL, SE, GR, RU.

Compatible with:

- Physiocontrol
- Lifepak 500 AED pink socket
- Lifepak 1000 AED
- Lifepak CR plus AED
- Lifepak Express AED