



Product Code	55020
Unit of sale	10 pairs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
--------	-----

RDM (NSIS)	2468982
CND	C020401
GTIN	28033003260313
GMDN	47055

## Description

ADULT Compatible pads - disposable

A range of compatible pads BSI certified for patients >8 years (25 kg) for most defibrillators on the market.

Excellent adherence and adequate gel thickness.

Includes a couple of pads, cable and connector.

Maximum shelf life 30 months (24 months for 33584) from production date.

Latex free.

Compliance to regulations: Directive 93/42/EEC and amendments. Class IIb medical device.

ANSI/AAMI: DF-80; IEC/EN: 60601-1, 60601-2-4; ISO: 10993-1.

Multilingual envelope: GB, FR, IT, ES, PT, DE, NL, SE, GR, RU.

Made in Italy.

Compatible with:

- Agilent-Philips
  - Heartstream/HeartStart XL M4735A
  - Heartstream/HeartStart XLT M3500B
  - HeartStart XL 861290
  - HeartStart 4000 M5500B
  - (all with Philips cable M3508A)
- Philips Laerdal Medical
  - Heartstream Fore Runner I AED models E, S, EM
  - Heartstart FR1 AED 940010XX - 940020XX - 940030XX
  - Heartstart FR2 AED M3841A - M3861A
  - Heartstart FR2+ AED, with ECG module M3840A, M3860A
  - Codemaster M1722A/B, M1723A/B, M1724A, M2475B
  - (with Philips cable M3507A and pad adapter 05-200)
  - Heartstart MRx M3535A (with Philips cable M3508A)