



| | |
|---------------|----------------|
| Product Code | 55022 |
| Unit of sale | 10 pairs |
| Minimum order | 1 |
| Type | Medical device |
| Class | II B |

| | |
|------------|---------|
| RDM (NSIS) | 2468987 |
| CND | C020401 |
| GMDN | 47055 |

Description

ADULT Compatible pads - disposable

A range of compatible pads BSI certified for patients >8 years (25 kg) for most defibrillators on the market.

Excellent adherence and adequate gel thickness.

Includes a couple of pads, cable and connector.

Maximum shelf life 30 months (24 months for 33584) from production date.

Latex free.

Compliance to regulations: Directive 93/42/EEC and amendments. Class IIb medical device.

ANSI/AAMI: DF-80; IEC/EN: 60601-1, 60601-2-4; ISO: 10993-1.

Multilingual envelope: GB, FR, IT, ES, PT, DE, NL, SE, GR, RU.

Made in Italy.

Compatible with:

- Agilent Hewlett-Packard
- Codemaster 100-XE-XL-XL+
- Heartstream XL, XLT, XL+

- Laerdal
- HeartStart 4000