



Product Code	55042
Unit of sale	10 pairs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
--------	-----

RDM (NSIS)	2469019
CND	C020401
GTIN	18054242447764
GMDN	47055

Description

ADULT Compatible pads - disposable

A range of compatible pads BSI certified for patients >8 years (25 kg) for most defibrillators on the market.

Excellent adherence and adequate gel thickness.

Includes a couple of pads, cable and connector.

Maximum shelf life 30 months (24 months for 33584) from production date.

Latex free.

Compliance to regulations: Directive 93/42/EEC and amendments. Class IIb medical device.

ANSI/AAMI: DF-80; IEC/EN: 60601-1, 60601-2-4; ISO: 10993-1.

Multilingual envelope: GB, FR, IT, ES, PT, DE, NL, SE, GR, RU.

Made in Italy.

Compatible with:

- GE
- Responder 1100 - 2000 - 3000 (direct connection)
- Cardioserv (direct connection)
- Dash Responder (direct connection)