

WOOD TONGUE DEPRESSOR - sterile (40 boxes of 50)



Product Code	25512
Unit of sale	box of 2000 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2815388
CND	V9001
EAN/UPC	8023279255126

Description

WOOD TONGUE DEPRESSOR - STERILE - box of 2,000 (40 boxes of 50 pcs)

Quality tongue depressors manufactured with selected high quality smooth birch wood. Edges accurately finished. Size: 150x18x1.6 mm. Individually packed in a paper envelope and sterilized.

Multilingual box: GB, FR, IT, ES, PT, DE, GR, Arabic