

PLASTIC TONGUE DEPRESSOR - adult/pediatric sterile
(20 box of 90)



Product Code	25516
Unit of sale	box of 1800 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	No
CH-REP	Yes

EAN/UPC	8023279255164
PARAF	987471099

Description

PLASTIC TONGUE DEPRESSORS ADULT/PEDIATRIC STERILE

- Tongue depressor - box of 1,800 (20 boxes of 90 pcs)

Translucent tongue depressors made of polystyrene with curved shape and smooth edges.
One end in adult size and one end in child size.
Latex and PVC free.

Curved shape allows perfect vision during examination avoiding that the hand of examiner obstruct the vision

- Child side 15 mm
- Adult side 20 mm
- Length 150 mm

Sterile model is individually packed in a peel-open pack (ethylene oxide).