

S/S STRAIGHT ARTERY FORCEPS - bubbles fantasy - 16 cm



Product Code	20552
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1630372
CND	L031306
EAN/UPC	802327920527
GTIN	0802327920527
GMDN	62468

Description

Stainless steel nurse straight artery forceps - 16 cm

• Model: Bubbles

Powder coated instruments are not autoclavable, can be sterilized in water at maximum 110°C