

# WIRE MESH BASKET for 35531-3



|               |                   |
|---------------|-------------------|
| Product Code  | 35538             |
| Unit of sale  | 1 pc              |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|         |               |
|---------|---------------|
| UK-REP  | Not required  |
| CH-REP  | Not required  |
| EAN/UPC | 8023279355383 |
| PARAF   | 903115956     |

## Description

Perforated Tray for Branson 8800 ultrasonic cleaner (codes 35531-35532-35533).