

GLASS BEAKER 600 ml



Product Code	35541
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279355413
PARAF	903115982

Description

Glass beaker 250 ml for Branson ultrasonic cleaner:

- Branson 2800 (codes 35501-35502-35503)
- Branson 3800 (codes 35510-35511-35512)
- Branson 5800 (codes 35520-35521-35522)
- Branson 8800 (codes 35531-35532-35533)