

GIMA BL3 WIDE SCREEN DIGITAL THERMOMETER °C - hang box



Product Code	25554
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1753078
CND	V03010102
EAN/UPC	8023279255546
UMDNS	14-032
PARAF	973922875

Description

BL3 WIDE SCREEN DIGITAL THERMOMETER - blister, flexible tip, water-proof, red/orange/green backlight after measurement, wide display 29.3 x 15.3 mm.

Low battery indicator. Battery 3 V button size (CR2032).

Display with three different illuminations:

- Green screen: temperature lower than 37.1 ° C
- Orange screen: temperature between 37.2 and 37.9 ° C
- Red screen (FEVER): temperature from 38.0 ° C up to 43.0 ° C

Provide accurate and quick read-out. Safe and easy to use with clear and large display.

- accuracy $\pm 0.1^{\circ}\text{C}$
- memory: last reading
- beeper function and auto shut-off
- plastic case
- multilanguage manual and box (GB, FR, IT, ES, PT, GR)