

TREVIRA CURTAINS for wing screen - white



Product Code	45572
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279455724

Description

CURTAINS IN TREVIRA 45xh 129 cm

Trevira curtains contain a permanent silver-based antimicrobial agent which is not affected by frequent washing.
UV resistant and flame-retardant material.
Disinfectable and washable at 60°.
Available in 6 colours.
Custom-made curtains available on request.

Multilingual manual: GB, FR, IT.

Made in Italy.

Made in Italy.