



|               |                |
|---------------|----------------|
| Product Code  | 25575          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | II A           |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | No  |

|            |               |
|------------|---------------|
| RDM (NSIS) | 894174        |
| CND        | V03010102     |
| EAN/UPC    | 8023279255751 |
| PARAF      | 977810720     |

Description

PROFESSIONAL INFRARED EAR THERMOMETER

- special design for hospital use
- durable probe for cover replacement up to 100,000 times
- probe cover ejector and probe cover detector
- one-second response time: fast and easy
- neck strap and 1 probe cover included
- fever alarm
- selectable °C/°F

Needs probe covers for use

Large LCD display

Multilingual box and manual: GB, FR, IT, ES.

Technical Specifications

Measurement range: 34~42.2 °C (93.2-108 °F)  
Accuracy: ±0.2 °C (0.4 °F) during 35.5~42 °C (95.9~107.6 °F), ±0.3 °C (0.5 °F) for other range  
Memory: 25 sets  
Buzzer: Yes  
Battery: CR2032 1 (Gima code 32212)  
Dimension - weight: 144.82x39.84x48.50 mm - 68.8 g