

VET RECTAL THERMOMETER



Product Code	25588
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	802327925581

Description

VET RECTAL THERMOMETER

Designed to be ergonomically correct for animals.
Setting for pets, pigs, cattle and other animals.

- Features:
- instant temperature in 6 seconds
 - beeps indicator
 - LCD display and memory
 - automatic shut-off after 15 seconds
 - display temperature in both Fahrenheit or Celsius scale
 - reliable, safe and easy to use

Multilingual box and manual: GB, FR, IT, ES, PT, GR.

Technical Specifications

Measurement range: 32~44°C (89.6-111.2 °F)
Accuracy: ± 0.1 °C (35~39°C)
Battery life is approx. 5,000 readouts
Dimension: 190x110x35 mm
Battery: 2x1.5V AAA
Measurement time: 6 sec.