



| | |
|---------------|----------------|
| Product Code | 55990 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| UK-REP | No |
| CH-REP | No |
| RDM (NSIS) | 2799008 |
| CND | V0899 |
| EAN/UPC | 6974580878059 |

Description

TRACTION SPLINT - adult

Traction splint relieves tissue pressure on fractured leg by avoiding further haemorrhages, shock and secondary damages to nerves, bones, muscles and blood vessels.

Adjustable length and automatic traction mechanism.

Made of aluminium alloy and nylon leg straps.

Supplied with carrying bag.

Size: 134x19x21 cm.

Folded size: 93x20x7 cm.

Weight: 2 kg.

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic