


GIMA STRETCHER 2 - foldable in 2 - black



Product Code	56005
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2799006
CND	V0899
EAN/UPC	6974580871210

Description

GIMA STRETCHER 2 - foldable in 2 - black

Foldable emergency stretcher made of high strength aluminium alloy.
Lightweight, small-sized, easy to carry and sterilize.

Supplied with 2 safety belts made of polypropylene fibre with plastic buckle.

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic

Technical Specifications

Open size cm: 185x50x3
Folded size cm: 92.5x50x7
Weight: 6.5 kg
Load: kg 159 kg