

MODULAR TROLLEY painted steel with 2+1+1 drawers



Product Code	45602
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	No
--------	----

CH-REP	No
--------	----

EAN/UPC	8023279456028
---------	---------------

GTIN	08023279456028
------	----------------

Description

MODULAR TROLLEYS SYSTEM

- Trolley with 2 small, 1 medium, 1 large drawers
- Painted steel

Modular trolleys made of painted steel or stainless steel, one-piece structure, with wide technopolymer plastic worktop with borders on 3 sides to be cleaned easily, no sharp corners, ergonomic tubular push-handle. Different drawers (small-yellow, medium-blue, large-red), depending on model, with epoxy painted panels and integrated handle grip, silent sliding with end-run block. Possibility to remove the drawers from the guides for an easy cleaning and disinfection. Central lock key. Overall bumpers and 4 antistatic Ø 125 mm castors, 2 with brakes.

Size: 67x63x104 cm.

Weight: 53 kg.

Made in Italy.