



Product Code	56030
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2799006
CND	V0899
EAN/UPC	6974580870411

Description

GIMA SCOOPING STRETCHER - aluminium - silver

Gima scooping stretchers allow operators to raise patient without having to move the body, minimising the risk of further injury. Detachable design for transferring and immobilize fractured and severely injured patients. The clutch hinges at both ends of the stretcher allow to split the main frame vertically into two parts.

- mainly used for hospitals, sport centres and ambulances
- adjustable length according to patient's body
- bottom end (feet) has a narrow frame structure

2 safety belts made of polypropylene fiber with plastic buckle.

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic

Technical Specifications

Open size cm: 168/210x44x7
 Folded size cm: 120x44x9
 Weight: 7.4 kg
 Max. load: 160 kg