



Product Code	56035
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2799006
CND	V0899
EAN/UPC	6974580874044

Description

GIMA SCOOPING STRETCHER - HPDE - yellow

Gima scooping stretchers allow operators to raise patient without having to move the body, minimising the risk of further injury. Detachable design for transferring and immobilize fractured and severely injured patients. The clutch hinges at both ends of the stretcher allow to split the main frame vertically into two parts.

- mainly used for hospitals, sport centres and ambulances
- adjustable length according to patient's body
- bottom end (feet) has a narrow frame structure
- thermal stability of polymer surface for better comfort

3 safety belts made of polypropylene fiber with plastic buckle.

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic

Technical Specifications

Open size cm: 161/210x45x6
 Folded size cm: 120x44.5x7
 Weight: 8 kg
 Max. load: 160 kg