

DELTA LED OTOSCOPE PENLIGHT - black



Product Code	25625
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279256253
PARAF	980418166

Description

DELTA "LED" OTOSCOPE PENLIGHT - black

Delta penlights are supplied in aluminium box with 2 AAA batteries and holder for tongue depressors, otoscope head, 3 ear funnels and 5 black plastic tongue depressors.

Box and instructions in: GB, FR, IT, ES, PT, DE.