



| | |
|---------------|-------------------|
| Product Code | 21563 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

| | |
|---------|---------------|
| EAN/UPC | 8023279215632 |
|---------|---------------|

Description

JACKET WITH STUDS 60% COTTON - 40% POLYESTER

- Size: M
- Colour: white
- For woman
- Five hidden press studs
- Short sleeves
- 2 front and 1 breast pockets

Professional and stylish line with all sizes for men and women, made of 60% high quality cotton and 40% polyester.
Stock in 5 colours, other 22 colours available in 60 days.
White jacket 200 g/m2, coloured jacket 185 g/m2.

Technical Specifications

Composition: 60% cotton - 40% polyester
Weight: 200 g/m2 (white), 185 g/m2 (coloured)
Yarn count: warp: Ne 20 ± 2 weft: Ne 18 ± 2
Construction: warp: 40 ± 3 weft: 19 ± 3
Breaking strength - warp: 100 kg ±3%
Breaking strength - weft: 39 kg ±3% (white) 49 kg ±3% (coloured)
Fastness to washing: 4
Treatments and finishing: crease resist
Dimensional stability - washing at 90°C: 3%