

HEINE MINI-C EXAMINATION LED LIGHT - D008.73.109



Product Code	25630
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
--------	-----

EAN/UPC	4053755198535
UMDNS	12-276
GMDN	36843

Description

HEINE MINI-C examination light LEDHQ - D008.73.109

Compact, long-life pocket diagnostic light.
20,000 switch cycles guaranteed.
Replaceable batteries size AAA (included).

Made in Germany.