

"ELEGANCE" TORCH - metal - blue



Product Code	25634
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
--------	----

EAN/UPC	8023279256345
PARAF	903120436
GMDN	12347

Description

Professional penlights, also suitable for diagnostic, with bulb or white LED with quick clip shifter mechanism or push button (25634). Argenta and Elegance made of anodized aluminium, Nova in topquality ABS plastic housing. Supplied with 2x1.5V AAA batteries.

Elegance torch
Light: bulb
Light intensity: 4.5 lux
Material: metal
Supplied with batteries: 2 x AAA
Colour: blue
ON/OFF by: Push button

Made in Germany.