

MASKOP MASK with wrap-around visor - ear loops



Product Code	25636
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2218687
CND	T020604
EAN/UPC	8030413033434
PARAF	987472406

Description

MASK-OP® with wrapping visor
- ear loops
- green

Liquid splash resistant and anti-fogging surgical mask with antistatic visor.

- 4 layers: three made of non-allergenic woven nonwoven and one of polyethylene net
- nuisance-free and adaptable nose-clip
- filtering efficiency (B.F.E.) =99%, according to EN 14683:2019+AC:2019
- resistance to breath <30,0 Pa/cm2, according to EN 14683:2019+AC:2019
- latex-free
- non sterile
- liquid splash resistance: resistant at pressure of 16 kPa, according to ISO 22609, related to EN 14683:2019+AC:2019
- anti-fogging internal stripe made of expanded polyethylene
- Bioburden = 30 CFU/g, according to EN 14683:2019+AC:2019

Medical Device Class I Regulation EU 2017/745.

Multilingual box: GB, FR, IT, ES, PT, DE, GR, PL, RU, Arabic.

Made in Italy.