

"ARGENTA" TORCH - metal - green



Product Code	25638
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
--------	----

RDM (NSIS)	2755589
CND	V9099
EAN/UPC	8023279256383
PARAF	903120451
GMDN	12347

Description

Professional penlights, also suitable for diagnostic, with bulb or white LED with quick clip shifter mechanism or push button (25634). Argenta and Elegance made of anodized aluminium, Nova in topquality ABS plastic housing. Supplied with 2x1.5V AAA batteries.

Argenta torch
Light: bulb
Light intensity: 4.5 lux
Material: metal
Supplied with batteries: 2 x AAA
Colour: green
ON/OFF by: Pocket clip

Made in Germany.