

4 PRACTI-MAN ADVANCE CPR MANIKINS



Product Code	56501
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8437010778162

Description

4 PRACTI-MAN ADVANCE CPR MANIKINS

CPR manikins for massive training events:

- adult and child CPR training with a single manikin
- 3-position selector dial: adult, child and test
- real-like system for breathing function
- jaw thrust maneuver
- AED training

Includes: 4 Practi-man, 20 lungs, 8 valves, 1 transport bag with mat.

Multilingual manual: GB, FR*, IT, ES*, PT*, DE*
(*Only on www.gimaitaly.com)