



Product Code	56509
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	7045432094862

Description

LAERDAL LITTLE BABY QCPR MANIKIN - 133-01050

Provided with QCPR feedback technology to help instructors to improve CPR training quality, efficiency and learner engagement. In compliance with current guidelines, they provide data about the rate and depth of chest compression, complete pressure relaxation, the time of CPR pause and ventilation, through the new QCPR Training Apps, skill guide available via RC cable.

It comes with:

- Little Baby QCPR Manikin
- 6 wipes
- manikin jumpsuit
- softpack carrying case
- user guide
- product information
- 2 AA batteries