

LAERDAL LITTLE ANNE QCPR - 120-01050



Product Code	56510
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	0704543207997

Description

LAERDAL LITTLE ANNE QCPR MANIKIN - 123-01050

Provided with QCPR feedback technology to help instructors to improve CPR training quality, efficiency and learner engagement. In compliance with current guidelines, they provide data about the rate and depth of chest compression, complete pressure relaxation, the time of CPR pause and ventilation, through the new QCPR Training Apps, skill guide available via RC cable.

- It comes with:
- Little Anne QCPR Manikin (manikin, 2 faces, 2 airways)
 - 6 manikin wipes
 - jacket
 - soft carrying case / training mat