



Product Code	56512
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	0704543209559

Description

LAERDAL LITTLE FAMILY PACK QCPR - 136-01050

Little Family Pack is a convenient package of agespecific CPR training manikins (Little Anne QCPR, Little Junior QCPR and Little Baby QCPR) with a wheeled carrying case for ease of transport and storage.

All are provided with QCPR feedback technology to help instructors to improve CPR training quality, efficiency and learner engagement.

In compliance with current guidelines, they provide data about the rate and depth of chest compression, complete pressure relaxation, the time of CPR pause and ventilation, through the new QCPR Training Apps, skill guide available via RC cable.

It comes complete with: Little Anne QCPR Manikin (manikin, 2 faces, 2 airways), Little Junior QCPR Manikin (manikin, 2 faces, 2 airways), Little Baby QCPR Manikin (manikin, 6 airways), 18 manikin wipes, 3 training mats and a large bag on wheels.