

TUBOLARE LATEX TOURNIQUET



Product Code	25733
Unit of sale	box of 200 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2249000
CND	V9003
EAN/UPC	8023599002042
PARAF	964240271

Description

GIMA TOURNIQUETS

- Model: Tubolare Tourniquet
- Material: latex
- Colour: yellow
- Length: 33 cm
- Fastening system: manual

Efficient and simple solutions for professional medical use.
Latex: Ø 5 mm, for manual use