

GASKET FOR H100 with serial number A,D,E,F,K,L

GUARNIZIONE / Gasket	
A D E F K L	☐
}= 8 45 1045 JUN 0SB	
B C	☐
}= 8 45 1045 099 0SB	
G	☐
}= 8 45 1045 DE3 0SB	
H	☐
}= 8 45 1045 MAG 0SB	
J	☐
}= 8 45 1045 M99 0SB	

ETICHETTA MATRICOLA / Matriculation label

Model

.....

Code

1 11 1105 000

Mtr.

99 ... A 0151

Year

1999

Max kPa

200

lt.

9

Max °C

130

Hz

50

V

230~

A

F8

W

1500

CE

Made in Italy

Product Code	35747
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	
UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279357479

Description

Gasket for Gima vertical autoclave H100 (cod. 35746-35656).

WARNING:
This gasket is only compatible for the autoclaves H100 with the following serial numbers:
A - D - E - F - K - L (see image, advise if you have different serial number).