

DISPOSABLE PVC LARYNGEAL MASK N 2.5



Product Code	57583
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2565401
CND	R01020199
EAN/UPC	6972091921912

Description

DISPOSABLE PVC LARYNGEAL MASKS

Laryngeal Mask N 2.5

Patient weight: 20 - 30 kg
Cuff inflation volume: 14 ml
Tube length: 130 mm
Minimum order: box of 50

Suitable for anesthesia and emergency medicine for airway management, causing less pain and coughing.

Tube with an inflatable cuff that is inserted into the pharynx.

Mask is made of PVC and consists of airway tube, laryngeal mask, connector, inflating tube, valve, pilot balloon, deflating flake and annectent back.

- allow rapid access
- provide airway for spontaneous or controlled ventilation
- tolerated at lighter anesthetic planes
- sterile and latex free.