



Product Code	57586
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2565401
CND	R01020199
EAN/UPC	6972091921943

## Description

### DISPOSABLE PVC LARYNGEAL MASKS

Laryngeal Mask N 5.0  
 Patient weight: 70 - 100 kg  
 Cuff inflation volume: 40 ml  
 Tube length: 195 mm  
 Minimum order: box of 50

Suitable for anesthesia and emergency medicine for airway management, causing less pain and coughing.

Tube with an inflatable cuff that is inserted into the pharynx.

Mask is made of PVC and consists of airway tube, laryngeal mask, connector, inflating tube, valve, pilot balloon, deflating flake and annectent back.

- allow rapid access
- provide airway for spontaneous or controlled ventilation
- tolerated at lighter anesthetic planes
- sterile and latex free.