

TROLLEY FOR PATIENT MONITOR - needs adaptor plate



Product Code	45760
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2329238
CND	V0880
EAN/UPC	8053736712985
UMDNS	10635
GMDN	36801

Description

TROLLEY FOR MONITORS - without plate adapter

Trolley specially designed for support and mobility of patient monitors and ventilators of most manufacturers.

Trolley needs plate adapter (not included) to be chosen depending on brand/model.

- 10 mm thickness steel base with high resistance ABS cover
- column in extruded aluminium with epoxy powders oven coated surface treatment
- ergonomic rear handle in ABS with reinforced structure
- 5 twin pivoting anti-noise castors Ø 65 mm made of thermoplastic rubber and ABS, 3 with brake

Storage basket made of ABS plastic (30.5x15.6xh 14.9 cm).

Size: 52.5x50xh 105 cm.

Max load: 15 kg.

Made in Italy.