

UNCUFFED ENDOTRACHEAL TUBE diam. 3 mm



Product Code	57640
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2567174
CND	R010399
EAN/UPC	6972091920199

Description

PVC UNCUFFED ENDOTRACHEAL TUBES - sterile

Uncuffed Endotracheal tubes Ø 3.0 mm

Tube length: 17.0 cm

Minimum order: box of 100

- bull-nose tips
- split-resistant radiopaque lines
- smooth Murphy eyes
- kink-resistant thermosensitive tubes