

UNCUFFED ENDOTRACHEAL TUBE diam. 4 mm



Product Code	57641
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2567174
CND	R010399
EAN/UPC	6972091920212

Description

PVC UNCUFFED ENDOTRACHEAL TUBES - sterile

Uncuffed Endotracheal tubes Ø 4.0 mm
Tube length: 20.5 cm
Minimum order: box of 100

- bull-nose tips
- split-resistant radiopaque lines
- smooth Murphy eyes
- kink-resistant thermosensitive tubes