

UNCUFFED ENDOTRACHEAL TUBE diam. 4.5 mm



Product Code	57642
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2567174
CND	R010399
EAN/UPC	6972091920229

Description

PVC UNCUFFED ENDOTRACHEAL TUBES - sterile

Uncuffed Endotracheal tubes Ø 4.5 mm

Tube length: 22.5 cm

Minimum order: box of 100

- bull-nose tips
- split-resistant radiopaque lines
- smooth Murphy eyes
- kink-resistant thermosensitive tubes