


# CUFFED ENDOTRACHEAL TUBE diameter 3 mm



Product Code	57650
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2565429
CND	R010399
EAN/UPC	6972091920021

## Description

### PVC CUFFED ENDOTRACHEAL TUBES - sterile

Cuffed Endotracheal tubes Ø 3.0 mm  
Tube length: 17.0 cm  
Minimum order: box of 100

Manufactured with strong tube walls to ensure the proper ventilation for all patients.  
A variety of sizes ensures a safe airway is maintained throughout the surgery.

- high-volume, low-pressure, barrel-shaped cuffs
- responsive pilot balloons
- finely heat-welded cuffs
- bull-nose tips
- split-resistant radiopaque lines
- smooth Murphy eyes
- kink-resistant thermosensitive tubes