


CUFFED ENDOTRACHEAL TUBE diameter 8.5 mm



Product Code	57661
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2567174
CND	R010399
EAN/UPC	6972091920137

Description

PVC CUFFED ENDOTRACHEAL TUBES - sterile

Cuffed Endotracheal tubes Ø 8.5 mm
Tube length: 32.5 cm
Minimum order: box of 100

Manufactured with strong tube walls to ensure the proper ventilation for all patients.
A variety of sizes ensures a safe airway is maintained throughout the surgery.

- high-volume, low-pressure, barrel-shaped cuffs
- responsive pilot balloons
- finely heat-welded cuffs
- bull-nose tips
- split-resistant radiopaque lines
- smooth Murphy eyes
- kink-resistant thermosensitive tubes