

CUFFED TRACHEOSTOMY TUBE diam. 7.5 mm



Product Code	57682
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2724821
CND	R0105020102
EAN/UPC	6972091921271

Description

Cuffed tracheostomy tubes Ø 7.5 mm

Curve tubes designed to rapidly provide an airway within seconds during an emergency by inserting them into a tracheostomy stoma, a hole created in the neck and windpipe.

The hollow tube whether with or without cuff, is directly placed into the trachea through a surgical incision, using a wire-guided progressive dilation technique.

- tube is made of medical grade, non-toxic, clear and kinkresistant PVC to protect delicate mucosal tissue
- smooth rounded obturator's tip reduces trauma during insertion
- neck plate with printed size for easy identification
- standard connector
- radio-opaque line helps to determine exact tube location

Tube includes: cuff, inflation line, valve, pilot balloon, obturator, neck plate, connector and tube ties.

Sterile, single patient. Class IIa. DEHP-free.

Made with no natural rubber or latex.

Cuffed tracheostomy tubes additional features

- high-volume low-pressure cuff provides an effective lowpressure seal
- soft, thin-walled cuff ensures effective sealing and atraumatic intubation and extubation
- efficient one-way valve for cuff inflation and deflation