

WASTE CONTAINER 50 l - single lid



Product Code	25776
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8013183125660
PARAF	977810795

Description

WASTE CONTAINERS - 3 sizes
Containers for disposal of special hospital waste in 3 sizes, with 2 different lids.

SC LINE WASTE CONTAINERS

- Tide: Container 50 l
 - Lid: single
 - Max load kg: 20
 - Size: 40x30xh 56 cm
- extremely solid and shock resistant at high and low temperature
 - hermetic seal by peripheral adhesive on the lid (single or double)
 - lateral and central handles tested for twice the max weight
 - stacking up to a maximum of 3 m with full load
 - conformity: ISO 23907-1:2019

Made in Italy.