

# WASTE CONTAINER 60 I - double lid



Product Code	25779
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8013183134471
PARAF	977810821

## Description

### WASTE CONTAINERS - 3 sizes

Containers for disposal of special hospital waste in 3 sizes, with 2 different lids.

#### SC LINE WASTE CONTAINERS

- Tide: Container 60 I
- Lid: double
- Max load kg: 25
- Size: 40x30xh 65 cm

- extremely solid and shock resistant at high and low temperature
- hermetic seal by peripheral adhesive on the lid (single or double)
- lateral and central handles tested for twice the max weight
- stacking up to a maximum of 3 m with full load
- conformity: ISO 23907-1:2019

Made in Italy.